


PRAMA

Product Overview



Neck morphology

Prama, the one and only intramucosal implant, was designed with a converging neck to maximize the thickness of soft tissues. To answer all clinical situations, Prama implants are available in three different heights of the intramucosal neck:



Short Neck
neck h 1.80 mm



Regular Neck
neck h 2.80 mm



Long Neck
neck h 3.80 mm

Prama Slim

Prama Slim was specifically developed for those areas with limited horizontal dimension. This implant has a straight intramucosal neck 1.80 mm or 2.80 mm high.



ø 3.30 Short Neck
neck h 1.80 mm



ø 3.30 Regular Neck
neck h 2.80 mm



Note: The use of Prama Slim for single crowns is just intended for the upper lateral incisors and for the lower central and lateral incisors. In all the other positions it can only be used for multiple rehabilitations.

Neck morphology and prosthesis

The concept of intramucosal Implant is very simple from the prosthetic point of view: considering the soft tissues thickness, it will be chosen the height of the neck between 1.80 mm, 2.80 mm and 3.80 mm and it will be decided whether to use a prosthesis that lays on the implant platform or whether to partially embrace the converging neck with dedicated Prama IN components.

All three necks can receive the prosthesis that closes at the platform level, while the Prama IN prosthesis that embraces the neck for 0.50 mm can only be used with 2.80 mm and 3.80 mm necks.

The Prama IN components that embrace the neck for 1.50 mm in height can only be used on the 2.80 mm neck.

		Standard prosthesis laying on the implant platform	0.50 mm Prama IN prosthesis	1.50 mm Prama IN prosthesis
Prama Slim				
	Neck			
	1.80  ø 3.30	✓	✗	✗
	2.80  ø 3.30	✓	✗	✗
Prama, Prama RF and Prama RF SL	1.80  ø 3.80 ø 4.25	✓	✗	✗
	2.80  ø 3.80 ø 4.25 ø 5.00	✓	✓	✓
	3.80  ø 3.80 ø 4.25 ø 5.00	✓	✓	✗

Guide to implant codes

Prama with cylindrical morphology

Morphology and neck	Surface	Diameter	Length
LAS: 1.80 mm - Short Neck LA: 2.80 mm - Regular Neck LAL: 3.80 mm - Long Neck	ZT: ZirTi surface and UTM neck	330: 3.30 mm - Prama Slim 380: 3.80 mm 425: 4.25 mm 500: 5.00 mm	060: 6.00 mm 085: 8.50 mm 100: 10.00 mm 115: 11.50 mm 130: 13.00 mm 150: 15.00 mm
example: LAS-ZT-380-100			
LAS-	ZT-	380-	100

Prama with conical morphology

Morphology and neck	Surface	Diameter	Spira		Length
			Regular	Large	
LSS: 1.80 mm - Short Neck LS: 2.80 mm - Regular Neck LSL: 3.80 mm - Long Neck	ZT: ZirTi surface and UTM neck	380: 3.80 mm 425: 4.25 mm 500: 5.00 mm	-	SL	060: 6.00 mm 085: 8.50 mm 100: 10.00 mm 115: 11.50 mm 130: 13.00 mm 150: 15.00 mm
example: LSL-ZT-380SL-100					
LSL-	ZT-	380		SL-	100

Prama Slim - Short and Regular Neck



h	ø 3.30
6.00	
8.50	LAS-ZT-330-085
10.00	LAS-ZT-330-100
11.50	LAS-ZT-330-115
13.00	LAS-ZT-330-130
15.00	LAS-ZT-330-150

Prama Slim
Short Neck

Surgical cover screw	
	L-VT-340



h	ø 3.30
6.00	LA-ZT-330-060
8.50	LA-ZT-330-085
10.00	LA-ZT-330-100
11.50	LA-ZT-330-115
13.00	LA-ZT-330-130
15.00	LA-ZT-330-150

Prama Slim
Regular Neck

Surgical cover screw	
	L-VT-340

Prama - Short Neck



h	ø 3.80	ø 4.25
6.00		
8.50	LAS-ZT-380-085	LAS-ZT-425-085
10.00	LAS-ZT-380-100	LAS-ZT-425-100
11.50	LAS-ZT-380-115	LAS-ZT-425-115
13.00	LAS-ZT-380-130	LAS-ZT-425-130
15.00	LAS-ZT-380-150	LAS-ZT-425-150

Prama

Surgical cover screw	
	L-VT-340



h	ø 3.80	ø 4.25
6.00		
8.50	LSS-ZT-380-085	LSS-ZT-425-085
10.00	LSS-ZT-380-100	LSS-ZT-425-100
11.50	LSS-ZT-380-115	LSS-ZT-425-115
13.00		
15.00		

Prama RF

Surgical cover screw	
	L-VT-340

Prama - Regular Neck



Prama



Prama RF



Prama RF SL

h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00	LA-ZT-380-060	LA-ZT-425-060	LA-ZT-500-060
8.50	LA-ZT-380-085	LA-ZT-425-085	LA-ZT-500-085
10.00	LA-ZT-380-100	LA-ZT-425-100	LA-ZT-500-100
11.50	LA-ZT-380-115	LA-ZT-425-115	LA-ZT-500-115
13.00	LA-ZT-380-130	LA-ZT-425-130	LA-ZT-500-130
15.00	LA-ZT-380-150	LA-ZT-425-150	LA-ZT-500-150

h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00		LS-ZT-425-060	LS-ZT-500-060
8.50	LS-ZT-380-085	LS-ZT-425-085	LS-ZT-500-085
10.00	LS-ZT-380-100	LS-ZT-425-100	LS-ZT-500-100
11.50	LS-ZT-380-115	LS-ZT-425-115	LS-ZT-500-115
13.00	LS-ZT-380-130	LS-ZT-425-130	LS-ZT-500-130
15.00	LS-ZT-380-150	LS-ZT-425-150	LS-ZT-500-150

h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00			
8.50	LS-ZT-380SL-085	LS-ZT-425SL-085	LS-ZT-500SL-085
10.00	LS-ZT-380SL-100	LS-ZT-425SL-100	LS-ZT-500SL-100
11.50	LS-ZT-380SL-115	LS-ZT-425SL-115	LS-ZT-500SL-115
13.00	LS-ZT-380SL-130	LS-ZT-425SL-130	LS-ZT-500SL-130
15.00	LS-ZT-380SL-150	LS-ZT-425SL-150	LS-ZT-500SL-150

[Surgical cover screw](#)



L-VT-340

Prama - Long Neck



Prama



Prama RF



Prama RF SL

h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00	LAL-ZT-380-060	LAL-ZT-425-060	LAL-ZT-500-060
8.50	LAL-ZT-380-085	LAL-ZT-425-085	LAL-ZT-500-085
10.00	LAL-ZT-380-100	LAL-ZT-425-100	LAL-ZT-500-100
11.50	LAL-ZT-380-115	LAL-ZT-425-115	LAL-ZT-500-115
13.00	LAL-ZT-380-130	LAL-ZT-425-130	
15.00			

h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00		LSL-ZT-425-060	LSL-ZT-500-060
8.50	LSL-ZT-380-085	LSL-ZT-425-085	LSL-ZT-500-085
10.00	LSL-ZT-380-100	LSL-ZT-425-100	LSL-ZT-500-100
11.50	LSL-ZT-380-115	LSL-ZT-425-115	LSL-ZT-500-115
13.00	LSL-ZT-380-130	LSL-ZT-425-130	
15.00			

h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00			
8.50	LSL-ZT-380SL-085	LSL-ZT-425SL-085	LSL-ZT-500SL-085
10.00	LSL-ZT-380SL-100	LSL-ZT-425SL-100	LSL-ZT-500SL-100
11.50	LSL-ZT-380SL-115	LSL-ZT-425SL-115	LSL-ZT-500SL-115
13.00			
15.00			

[Surgical cover screw](#)



L-VT-340

Surgical kits

Prama, Prama RF and Prama RF SL surgical kit



ZPRAMA-INT complete surgical grommetless kit of the instruments necessary for Prama and Prama RF implants

L-TRAYL-INT radel instrument grommetless tray

Drilling kit for Prama cylindrical short implants




ZSHORTY-INT drilling kit Shorty for Prama short implants

SHORTY-TRAY-INT radel instrument tray

GROMMET-CA-1 spare grommets for drills or contra-angle instruments (5 units package)

Please note: the ZSHORTY-INT kit allows you to insert the Shorty implants of the Premium, Outlink² and Prama implant lines: please refer to the MC-IMP-PRAMA surgical manual for the correct sequence of use of the only instruments necessary for Prama diameters. ZSHORTY-INT kit does not allow the insertion of Prama RF short implants, for which the instruments contained in the L-INTEGRA-060 set at page 8 are necessary.

Surgical kits for guided insertion of Prama implants




ZGS-PRO-A-INT grommetless surgical kit including all the instruments for guided insertion of Prama implants

GS-PROA-INTEGRA-060 set of drills for the preparation of Prama implants, height 6.00 mm

GSPROA-TRAY-INT radel instrument grommetless tray

For further information on guided insertion of implants, please consult the surgical manual MC-UD-EPLANPROA-E.

Surgical kits for guided insertion of Prama RF and Prama RF SL implants



ZGS-PRO-S-INT grommetless surgical kit including all the instruments for guided insertion of Prama RF and Prama RF SL implants

GS-PROS-INTEGRA-060 set of drills for the preparation of Prama RF implants, height 6.00 mm

GSPROS-TRAY-INT radel instrument grommetless tray

For further information on guided insertion of implants, please consult the surgical manual MC-UD-EPLANPROS-E.

Instruments included in Prama surgical kit

Precision drill			
	FS-230		
Pilot drill and related stops			
	FPT3-200-LXS		
	h		
	8.50	STOP4-200-085	
	10.00	STOP4-200-100	
	11.50	STOP4-200-115	
	13.00	STOP4-200-130	
	15.00	STOP4-200-150	
Intermediate drills and related stops			
	FFT3-250-LXS	2.50 cylindrical drill	
	SE-FK250	2.50 conical drill	
	h	cylindrical	conical
	8.50	STOP4-250-085	CSR-STOP-2028-085
	10.00	STOP4-250-100	CSR-STOP-2028-100
	11.50	STOP4-250-115	CSR-STOP-2028-115
	13.00	STOP4-250-130	CSR-STOP-2028-130
	15.00	STOP4-250-150	CSR-STOP-2028-150
Prama RF Reply			
	h	ø 3.80 	ø 4.25  ø 5.00 
	8.50	SH-380-085-RP	SH-425-085-RP SH-500-085-RP
	10.00	SH-380-100-RP	SH-425-100-RP SH-500-100-RP
	11.50	SH-380-115-RP	SH-425-115-RP SH-500-115-RP
	13.00	SH-380-130-RP	SH-425-130-RP SH-500-130-RP
	15.00	SH-380-150-RP	SH-425-150-RP SH-500-150-RP
Prama bone taps with hexagonal connector			
	ø 3.30 	ø 3.80 	ø 4.25  ø 5.00 
	A-MS-330	A-MS-380	A-MS-425 A-MS-500
Prama RF bone taps with right angle shank			
	ø 3.80 	ø 4.25 	ø 5.00 
	SH-MS-380-CA	SH-MS-425-CA	SH-MS-500-CA
Screwdrivers			
	HSMX-20-DG	digital, extra short	
	HSM-20-DG	digital, short	
	HSML-20-DG	digital, long	
	HSM-20-EX	with hexagonal connector	
	HSML-20-EX	with hexagonal connector, long	
	HSM-20-CA	with right angle shank	
Dynamometric ratchet			
	CRI5-KIT		
Pins for neck positioning			
	h		
	8.50	L-PP-085	
	10.00	L-PP-100	
	11.50	L-PP-115	
	13.00	L-PP-130	
	15.00	L-PP-150	
Parallelism pin with depth laser marking			
	PPTL-2-28		

Prama final drills and related stops					
		ø 3.30 	ø 3.80 	ø 4.25 	ø 5.00 

Optional sets

L-INTEGRA-060 set

L-INTEGRA-060 set includes all the necessary components for the insertion of Prama and Prama RF h 6.00 mm.

Slots for these instruments are provided in the surgical kit.

Initial drills for Prama h 6.00 mm and related stop		
	FPS-200	pilot drill
	STOPS4-200-060	stop for pilot drill
	FGS-200/300	intermediate drill 200/300
	FGS-340/425	intermediate drill 340/425

Instruments included in the ZSHORTY-INT kit.

Final drills for Prama RF h 6.00 mm and related stops		
	ø 4.25 	ø 5.00 
	SE-FK410	SE-FK500
	SE-STOP-FK410-060	SE-STOP -FK500-060

Final drills for Prama h 6.00 mm and related stops			
	ø 3.80 	ø 4.25 	ø 5.00 
	FFS-300	FFS-340	FFS-425
	STOPS4-300-060	STOPS4-340-060	STOPS4-425-060

Instruments included in the ZSHORTY-INT kit.

Pins for neck positioning of Prama h 6.00 mm	
	L-PP-060

Reply for Prama RF h 6.00 mm	
	SE-410-060-RP
	SE-500-060-RP

L-INTEGRA set

L-INTEGRA set includes all the necessary components for the insertion of Prama Slim ø 3.30 mm, the countersink drills and the pins for neck positioning.

These components are included in the ultimate version of the Prama surgical kit, but not in the previous version of the surgical kit, thus they are still available for those who have the previous version.

Pins for neck positioning		
	h	
	8.50	L-PP-085
	10.00	L-PP-100
	11.50	L-PP-115
	13.00	L-PP-130
	15.00	L-PP-150

Countersink drills			
	ø 3.80 	ø 4.25 	ø 5.00 
	A-FCC-380	A-FCC-425	A-FCC-500

Instruments for Prama Slim implants		
	FFT3-280-LXS	final drill
	h	
	8.50	STOP4-280-085
	10.00	STOP4-280-100
	11.50	STOP4-280-115
	13.00	STOP4-280-130
	15.00	STOP4-280-150
	A-MS-330	bone tap with hexagonal connector

Optional instruments (not included in Prama surgical kit)

Intermediate drills

	FG-200/280XS	intermediate drill 200/280
	FG-330/425XS	intermediate drill 330/425

Drills for distal sectors

ø		ø	
2.00	FPT5-200-LXS	3.30	FFT5-330-LXS
2.80	FFT5-280-LXS	3.40	FFT5-340-LXS
2.90	FFT5-290-LXS	3.60	FFT5-360-LXS
3.00	FFT5-300-LXS	4.25	FFT5-425-LXS
3.20	FFT5-320-LXS	4.55	FFT5-445-LXS

Prama bone taps

with right angle shank			
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●
	A-MS-330-CA	A-MS-380-CA	A-MS-425-CA
with hexagonal connector, short			
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●
	A-MS-330	A-MS-380	A-MS-425

Prama RF bone taps with hexagonal connector

	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	SH-MS-380	SH-MS-425	SH-MS-500


Optional prosthetic screwdrivers

	HSMXL-20-EX	with hexagonal connector, extra long
	HSMXS-20-CA	with right angle shank
	L-HSM-EX	Full Head, with hexagonal connector, short
	L-HSML-EX	Full Head, with hexagonal connector, long
	L-HSMXL-EX	Full Head, with hexagonal connector, extra long
	L-HSM-CA	Full Head, with right angle shank
	AVV-CAP-TIT-1	instrument for maintaining of the titanium cap for ball attachments
	BASCC-EX	for ball attachments, with hexagonal connector
	AVV2-ABUT	for straight P.A.D. abutments
	8926-SW	for Locator abutment, short
	8927-SW	for Locator abutment, long
	DSPDCLH-24	for Dynamic Abutment, length 24 mm
	DSPDCLH-32	for Dynamic Abutment, length 32 mm

Mounter

	MOU-EX230	mounter
	CM2	mounter stop key

Prama osteotomes

	flat tip	
	E-OS-020-PP	tip ø 0.20
	E-OS-090-PP	tip ø 0.90
	concave tip	
	E-OS-160-PC	tip ø 1.60
	E-OS-200-PC	tip ø 2.00
	E-OS-240-PC	tip ø 2.40

Prama RF osteotomes

	flat tip	round tip
h	ø 3.80 ●	
8.50 - 10.00	SH-OS-380-100-PP	SH-OS-380-100-PR
11.50	SH-OS-380-115-PP	SH-OS-380-115-PR
13.00	SH-OS-380-130-PP	SH-OS-380-130-PR
15.00	SH-OS-380-150-PP	SH-OS-380-150-PR
	ø 4.25 ●	
8.50 - 10.00	SH-OS-425-100-PP	SH-OS-425-100-PR
11.50	SH-OS-425-115-PP	SH-OS-425-115-PR
13.00	SH-OS-425-130-PP	SH-OS-425-130-PR
15.00	SH-OS-425-150-PP	SH-OS-425-150-PR
	ø 5.00 ●	
8.50 - 10.00	SH-OS-500-100-PP	SH-OS-500-100-PR
11.50	SH-OS-500-115-PP	SH-OS-500-115-PR
13.00	SH-OS-500-130-PP	SH-OS-500-130-PR
15.00	SH-OS-500-150-PP	SH-OS-500-150-PR

Extraction tools for intraoral removing of implants

	BC-EX230	short
	BL-EX230	long

Depth gauge

	PROF3
--	-------

P.A.D. carrier

	PAD-CAR
--	---------

Torque wrench with control lever

	TWL
--	-----

Parallelism pins

	PPTS-2-28	with depth laser marking, small
	PP-2/28	without depth laser marking, standard
	PPS-2/3	for short implants

Spare O-ring

	ORING180-088
---	--------------

Tissues conditioning

Healing abutments laying on the implant platform

Straight healing abutments		
	transgingival height	unique platform
	2.00	A-TMG-330-2
	3.00	A-TMG-330-3
	5.00	A-TMG-330-5
	7.00	A-TMG-330-7
Healing abutments with coronal enlargement		
	transgingival height	unique platform
	2.00	A-TMTCR-330-2
	3.00	A-TMTCR-330-3
	transgingival height	unique platform
	2.00	A-TMGB-330-2
	3.00	A-TMGB-330-3
	4.00	A-TMGB-330-4

Temporary posts

Temporary posts		
	unique platform	
	A-PPF-330-EX	in REEF resin, repositionable
	A-PPF-330	in REEF resin, rotating
	A-MPSC-330	Simple in PEEK, with Gr. 5 titanium base, repositionable
	A-MPSCI-330-EX	Simple in Gr. 5 titanium base, repositionable
	A-MPSCI-330	Simple in Gr. 5 titanium base, rotating
	A-CTI-330	in titanium with castable sleeve, rotating
	A-CCI-S	spare castable sleeve for luting technique

VM2-180 fixation screw included.

Prima IN healing abutments (closing on the implant neck)

Prima IN healing abutments in titanium						
	transgingival height	closure on the neck	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	2.50	0.50*	L-TMG-380-25-05	L-TMG-425-25-05	L-TMG-500-25-05	
	3.50	0.50*	L-TMG-380-35-05	L-TMG-425-35-05	L-TMG-500-35-05	
	2.50	1.50**	L-TMG-380-25-15	L-TMG-425-25-15	L-TMG-500-25-15	
	3.50	1.50**	L-TMG-380-35-15	L-TMG-425-35-15	L-TMG-500-35-15	
Prima IN healing abutments in REEF resin						
	transgingival height	closure on the neck	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	6.00	0.50*	L-TMGPF-380-05	L-TMGPF-425-05	L-TMGPF-500-05	
	6.00	1.50**	L-TMGPF-380-15	L-TMGPF-425-15	L-TMGPF-500-15	

Healing abutments in REEF resin: VM2-180 fixation screw included.

*Only for Regular and Long Neck. | **Only for Regular Neck.

Impression and model phase

Pick-up transfers				
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	A-TRA-330	L-TRA-380	L-TRA-425	L-TRA-500

VTTRA2-180-15 screw included.

Pull-up transfers	
	unique platform
	A-TRAP-330

Analog							
Short Neck							
	ø 3.30 		ø 3.80 		ø 4.25 		
A-ANA-330		LS-ANA-380		LS-ANA-425			
Regular Neck							
	ø 3.30 		ø 3.80 		ø 4.25 		ø 5.00
A-ANA-330		L-ANA-380		L-ANA-425		L-ANA-500	
Long Neck							
			ø 3.80 		ø 4.25 		ø 5.00
		LL-ANA-380		LL-ANA-425		LL-ANA-500	

Posts

Premade posts

Straight premade posts			
	h	unique platform	
	6.00	L-MD-340-6	straight, repositionable
	8.00	L-MD-340-8	

Straight premade posts shoulderless			
	h	unique platform	
	15.00	A-MD-330-15	repositionable

Preparable posts

Standard preparable posts (laying on the implant platform)		
	unique platform	
	A-MF-330	straight emergence
		
	A-MFS-330	wide emergence

VM2-180 fixation screw included.

Prima IN preparable posts (closing on the implant neck)				
	closure on the neck	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	0.50*	L-MF-380-05	L-MF-425-05	L-MF-500-05
	1.50**	L-MF-380-15	L-MF-425-15	L-MF-500-15

L-VMS-180 fixation screw included. | * Only for Regular and Long Neck | ** Only for Regular Neck.

Castable posts with a metal base and entirely castable posts

Castable abutments with a metal base		
	unique platform	
	A-UCR-330	gold alloy base, rotating
	A-UCRCO-330	cobalt chrome base, rotating
	A-UCRCO-330-EX	cobalt chrome base, repositionable
	A-UCTR-330-EX	titanium base, repositionable
	A-UCR-330-EX	gold alloy base, repositionable
	A-CCUCR-330	spare castable PMMA sleeve, fixation screw not included

VM2-180 fixation screw included.

Prima IN castable posts with a metal base (closing on the implant neck)					
	closure on the neck	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	0.50*	L-UCRCO-380-05	L-UCRCO-425-05	L-UCRCO-500-05	cobalt chrome base, repositionable
	1.50**	L-UCRCO-380-15	L-UCRCO-425-15	L-UCRCO-500-15	cobalt chrome base, repositionable
	0.50*	L-UCRCO-380-05-ROT	L-UCRCO-425-05-ROT	L-UCRCO-500-05-ROT	cobalt chrome base, rotating
	1.50**	L-UCRCO-380-15-ROT	L-UCRCO-425-15-ROT	L-UCRCO-500-15-ROT	cobalt chrome base, rotating
	0.50*	L-UCR-380-05-ROT	L-UCR-425-05-ROT	L-UCR-500-05-ROT	gold alloy base, rotating
					
		A-CCUCR-380	A-CCUCR-425	A-CCUCR-500	spare castable PMMA sleeve, fixation screw not included

* Only for Regular and Long Neck | ** Only for Regular Neck. | L-VMS-180 fixation screw included.

For further information please consult the prosthetic manual MP-IMP-PRAMA-E.

Angled premade posts			
	transgingival height	unique platform	
	0.00	L-MA15-340*	angled at 15°, repositionable
	1.00	L-MA15-340-1	
	2.00	L-MA15-340-2	angled at 25°, repositionable
	1.00	L-MA25-340-1	
	2.00	L-MA25-340-2	

L-VMS-180 fixation screw included. | *L-MA15-340 posts: L-VM-180 fixation screw included.

Posts for luting technique

Straight premade posts			
	h pilastro	unique platform	
	6.00	L-MD-340-6-ROT	conical, rotating
	8.00	L-MD-340-8-ROT	

Vite L-VMS-180 inclusa.

Posts in titanium and in cobalt chrome - rotating		
	unique platform	
	L-CT-340-ROT	Gr. 5 titanium posts, h 14.00 mm, straight emergence
	L-CCRCO-340-ROT	Cobalt chrome titanium posts, h 14.00 mm, straight emergence
	L-CTR-340-ROT	Gr. 5 titanium posts, h 14.00 mm, anatomical emergence
	L-CCRCOR-340-ROT	Cobalt chrome posts, h 14.00 mm, anatomical emergence

L-VMS-180 fixation screw included.

Dynamic Abutment for angled screw hole		
	unique platform	
	PD3PKH330/CC	cobalt chrome base, repositionable
	PD3PKR330/CC	cobalt chrome base, rotating
	PD3PKH330/P	entirely castable, repositionable
	PD3PKR330/P	entirely castable, rotating

TPDH18L66 fixation screw not included.

Entirely castable posts		
	unique platform	
	A-CC-330-EX	repositionable
	A-CC-330	rotating

VM2-180 fixation screw included.

Abutment

P.A.D. abutment (Disparallel Screwed Prosthesis)

Straight P.A.D. abutments		
	transgingival height	unique platform
	1.50	A-PAD-AD330-15
	3.00	A-PAD-AD330-30
	4.00	A-PAD-AD330-40

P.A.D. components for overstructures

Protection caps		
	unique platform	
	PAD-CG	in Gr. 5 titanium
	PAD-CGP	in PEEK

PAD-VP-140 fixation screw included.

Impression and model phase		
	unique platform	
	PAD-CAP	POM cap for impression taking, rotating
	PAD-CAP-EX	POM cap for impression taking, repositionable
	PAD-TRA	Gr. 5 titanium Pick-up transfer, rotating
	PAD-TRA-EX	Gr. 5 titanium Pick-up transfer, repositionable
	PAD-ANA	Gr. 5 titanium analog for P.A.D. abutment

Pick-up transfer: PAD-VTRAL-140 fixation screw included.

P.A.D. prostheses for D.P.F. technique		
	unique platform	
	PAD-CT-LV	spare titanium sleeve
	PAD-CC-LV	spare castable centring device
	PAD-TR-LV	spare anti-escape plug
	PAD-ORING-LV	spare protection O-ring
	PAD-LV	complete pack of the prosthetic components for the "D.P.F." technique
	BARC	castable bar, l 5 cm, ø 2.20 mm

Titanium sleeves: PAD-VP-140 fixation screw included.

Abutment P.A.D. angolati		
	transgingival height	unique platform
	3.00	A-PAD-AA330-173
	5.00	A-PAD-AA330-175
	3.00	A-PAD-AA330-303
	5.00	A-PAD-AA330-305

PAD-VM-180 fixation screw included.

P.A.D. components for prosthesis relining and luting technique		
	unique platform	
	PAD-CP	PEEK sleeve, rotating
	PAD-CP-EX	PEEK sleeve, repositionable
	PAD-CT	Gr. 5 titanium sleeve, rotating
	PAD-CT-EX	Gr. 5 titanium sleeve, repositionable
	PAD-CCEM	PMMA castable sleeve for Gr. 5 titanium sleeves, fixation screw not included
	PAD-CC	PMMA castable sleeve, rotating
	PAD-CC-EX	PMMA castable sleeve, repositionable

Sleeves in PEEK, Gr. 5 titanium and PMMA: PAD-VP-140 fixation screw included.

Castable abutments with a metal base		
	unique platform	
	PAD-UC	PMMA castable sleeve, gold alloy base, rotating
	PAD-UCRCO	PMMA castable sleeve, cobalt chrome base, rotating

PAD-VP-140 fixation screw included.

PLAIN abutments

PLAIN abutments		
	transgingival height	unique platform
	2.00	A-PLAIN-ABU330-2
	3.00	A-PLAIN-ABU330-3
	4.00	A-PLAIN-ABU330-4

Components for PLAIN abutments		
	unique platform	
	A-PLAIN-CG330	healing cap
	A-PLAIN-CC330	castable sleeve
	A-PLAIN-CT330	titanium sleeve
	A-PLAIN-TRA-330	Pick-up transfer
	A-PLAIN-ANA-330	analog

Castable and titanium sleeves: A-PLAIN-VP200 fixation screw included.

Pick-up transfer: A-PLAIN-VTRA200 fixation screw included.

For further information please consult the prosthetic manual MP-IMP-PRAMA-E.

CAD-CAM prosthetic solutions

Scanbodies for laboratory and intraoral use

Implant scanbody		
unique platform		
	L-INT-CAMTRA-L	h 10.00 mm, in titanium

Implant scanbody: VM2-180 fixation screw included.

Abutment scanbody		
unique platform		
	PAD-INT-CAMTRA-L	for P.A.D. abutment, in titanium
	A-INT-CAMTRA-PLAIN-L	for PLAIN abutment, in titanium

P.A.D. abutment scanbody: PAD-VP-140 fixation screw included.
PLAIN abutment scanbody: A-PLAIN-VP200 fixation screw included.

Intermediate elements for luting

Standard T-Connect (laying on the implant platform)			
	cone height	unique platform	
	4.00	A-BASTZR-S-330-4	repositionable
	6.00	A-BASTZR-S-330-6	
	4.00	A-BASTZR-M-330-4	rotating
	6.00	A-BASTZR-M-330-6	

VM2-180 fixation screw included.

T-Connect con diverse altezze transmucose			
	transgingival height	cone height	unique platform
	1.00	4.00	A-TCR-S-330-4-1
	2.00	4.00	A-TCR-S-330-4-2
	3.00	4.00	A-TCR-S-330-4-3
	1.00	4.00	A-TCR-M-330-4-1
	2.00	4.00	A-TCR-M-330-4-2
	3.00	4.00	A-TCR-M-330-4-3

VM2-180 fixation screw included.

Prama IN T-Connect supports (closing on the implant neck at 0.50 mm)*					
	cone height	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	4.00	L-BASTZR-S-380-4-05	L-BASTZR-S-425-4-05	L-BASTZR-S-500-4-05	repositionable
	6.00	L-BASTZR-S-380-6-05	L-BASTZR-S-425-6-05	L-BASTZR-S-500-6-05	
	4.00	L-BASTZR-M-380-4-05	L-BASTZR-M-425-4-05	L-BASTZR-M-500-4-05	rotating
	6.00	L-BASTZR-M-380-6-05	L-BASTZR-M-425-6-05	L-BASTZR-M-500-6-05	

*Only for Regular and Long Neck. | L-VMS-180 fixation screw included.

T-Connect for single crown (closing on the implant neck at 0.50 mm)*					
	cone height	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	5.00	L-TCIN-380-05	L-TCIN-425-05	L-TCIN-500-05	repositionable
	5.00	L-TCIN-380-05-ROT	L-TCIN-425-05-ROT	L-TCIN-500-05-ROT	rotating

*Only for Regular and Long Neck. | For these components use the standard VM2-180 fixation screw or the angled screw hole TPDH18L66. Screws are not included.

T-Connect for single crown (closing on the implant neck at 1.00 mm)*					
	cone height	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	5.00	L-TCIN-380-10	L-TCIN-425-10	L-TCIN-500-10	repositionable

*Only for Long Neck. | For these components use the standard VM2-180 fixation screw or the angled screw hole TPDH18L66. Screws are not included.

Prama IN T-Connect supports (closing on the implant neck at 1.00 mm)*					
	cone height	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	5.00	LL-TCIN-380-10	LL-TCIN-425-10	LL-TCIN-500-10	repositionable

*Only for Long Neck. | For these components use the standard VM2-180 fixation screw or the angled screw hole TPDH18L66. Screws are not included.

Interfase Dynamica with support base and palatal screw hole			
	cone height	unique platform	
	4.00	IND3PKH330/TIA	repositionable
	4.00	IND3PKR330/TIA	rotating

TPDH18L66 fixation screw not included.

Interfase without support base and with palatal screw hole			
	cone height	unique platform	
	6.00	L-MDT-340-6	repositionable
	6.00	L-MDT-340-6-ROT	rotating
	10.00	L-MDT-340-10	repositionable
	10.00	L-MDT-340-10-ROT	rotating

For these components use the standard VM2-180 fixation screw or the angled screw hole TPDH18L66. Screws are not included.

Sleeves for wax-up modeling



cone height	ø 3.30	ø 3.80	ø 4.25	ø 5.00
4.00	A-CCBAS-330-4	A-CCBAS-380-4	A-CCBAS-425-4	A-CCBAS-500-4
6.00	A-CCBAS-330-6	A-CCBAS-380-6	A-CCBAS-425-6	A-CCBAS-500-6

Titanium sleeve for P.A.D. abutment



PAD-BASTZR

PAD-VP-140 fixation screw included.

Chairside T-Connect supports with cementing cone h 5.00 mm



transgingival height	unique platform
0.50	A-TCS-330S-C
1.50	A-TCS-330S-C-15
2.50	A-TCS-330S-C-25

VM2-180 fixation screw included.

Chairside Prama In T-Connect supports (closing on the implant neck at 0.50 mm)*



ø 3.80	ø 4.25	ø 5.00
L-TCS-380S-05	L-TCS-425L-05	L-TCS-500L-05

*Only for Regular and Long Neck. | L-VMS-180 fixation screw included.

Premilled posts

Titanium Prama IN premilled posts (closing on the implant neck at 0.50 mm)*



ø 3.80	ø 4.25	ø 5.00
L-CAMT-38005-12-SW	L-CAMT-42505-12-SW	L-CAMT-50005-12-SW

VM2-180 fixation screw included. | *Only for Regular and Long Neck.

Titanium premilled post (laying on the implant platform)



unique platform
A-CAMT-330-12-SW

VM2-180 fixation screw included.

Conico conometric prosthesis

Straight Conico posts



transgingival height	unique platform
0.50	A-MD-TS-330-05
1.00	A-MD-TS-330-10
2.00	A-MD-TS-330-20
3.00	A-MD-TS-330-30

VM2-180 fixation screw included.

Gr. 2 titanium bars, for prostheses on conometric posts*



DW-BARRA1.2 circular profile
l 150 mm, ø 1.20 mm



DW-BARRA1.5 circular profile
l 150 mm, ø 1.50 mm



DW-BARRA1.8 circular profile
l 150 mm, ø 1.80 mm



DW-BARRA1x3 rectangular profile
l 100 mm, 3x1 mm



DW-BARRA2x4 rectangular profile
l 100 mm, 4x2 mm

* 5 pcs packaging.

Angled Conico posts



transgingival height	unique platform	
2.00	A-MA07-TS-330-2	angolati a 7.5°



2.00	A-MA15-TS-330-2	angolati a 15°
------	-----------------	----------------



2.00	A-MA20-TS-330-2	angolati a 22.5°
------	-----------------	------------------

VM2-180 fixation screw included.

Conometric components



CAP2-TS-DEF cap for partial or full-arch fixed prosthesis



CAP2-TS-IND non-rotating cap for single element



CAP2-TS-REM cap for the removal



CAP2-TS-IMP PEEK transfer cap for direct impression on Conico abutment



ANA2-MD-TS Conico post analog



CAP2-TS-PIN cap with pin, to be used for intraoral welding

Accessories for Conico system



PAR-PP parallelometer



PAR-INS-AS parallelometer insert

Remember to always purchase the parallelometer with the insert relating to the implant connection on which it is to be used.

Clinical and laboratory prosthetic screws

Standard fixation screws

	VM2-180	
	VM2-180-10	10 pcs packaging

Conical support fixation screws

	L-VMS-180	
	L-VMS-180-10	10 pcs packaging

Full Head technology fixation screws

	L-VM-180	
	L-VM-180-10	10 pcs packaging

Fixation screws for Dynamic Abutment and Interfase Dynamica

	TPDH18L66	
---	-----------	--

Screws for Pick-up transfer

	VTRA2-180-15	h 15.00
	VTRA2-180-20	h 20.00
	VTRA2-180-MAN	for manual screwing

Recommended torque for tightening prosthetic components

8-10 Ncm	surgical cover screws, healing abutments
20-25 Ncm	prosthetic screws
25-30 Ncm	direct screwing on implants prosthetic components
8-10 Ncm	transfer screws

Screws for P.A.D. abutments

	PAD-VM-180	
	PAD-VM-180-10	10 pcs packaging
	PAD-VTRAL-140	for Pick-up transfer, long
	PAD-VTRA-140	for Pick-up transfer, standard
	PAD-VTRAL-140-MAN	for Pick-up transfer, manual screwing
	PAD-VP-140	spare screw
	PAD-VP-140-10	10 pcs packaging

Screws for PLAIN abutments

	A-PLAIN-VP200	for PLAIN abutments
	A-PLAIN-VP200-10	10 pcs packaging
	A-PLAIN-VTRA200	for PLAIN abutment transfers

Screws for CAD-CAM custom made prosthesis

	A-CAMTVABU180	for zirconium oxide single crowns and ECHO overstructures
	PAD-VCAM140	for zirconium oxide ECHO overstructures on P.A.D. abutments
	CAMPRON205-10	spare ring washers for the head of the screw

Warning

Fixation screws are always included in the packaging of each prosthetic component, unless otherwise indicated.
Please use laboratory screws for laboratory phases and keep the screw included with the prosthetic component for the final fastening in the patient's mouth.

Solutions for overdenture

Locator and Locator R-Tx abutments

Locator Abutment		
	transgingival height	
	1.00	1670
	2.00	1671
	3.00	1672
	4.00	-
	5.00	1674
	6.00	-

Accessories and spare parts for Locator abutments		
	unique platform	
	8515	low retention capacity for impression taking, black, in polyethylene, 4 pcs packaging
	8524	retention 5 lb = 2268 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	8527	retention 3 lb = 1361 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	8529	retention 1,5 lb = 680 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	8547	retention 4 lb = 1814 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8915	retention 2 lb = 907 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8548	retention 1 lb = 450 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8514	spacer rings in silicon rubber, 20 pcs packaging
	8519-2	kit: 2 Gr. 5 titanium caps, 2 spacer rings and 2 nylon retainers for each of the 4 different retention capacities
	8540-2	kit: 2 Gr. 5 titanium caps, 2 spacer and 2 nylon retainers for each of the 4 different retention capacities meant for disparallelisms
	8550-2	kit: 2 Gr. 5 steel caps, 2 spacer rings and 2 nylon retainers for each of the 4 different retention capacities

Accessories and spare parts for Locator abutments		
	unique platform	
	8530	4 aluminium analogs
	8505	4 aluminium transfers, 4 retainers cod. 8515 included
	8517	4 parallelism pins in nylon
	9530	steel plate AISI 316 L

Locator Core Tool		
	8393	complete Locator Core Tool
	8397	steel push rod for inserting the retainers in the caps
	8390	steel driver for abutment screwing/unscrewing
	8394	retention jacket for the driver, to transfer the Locator abutments into the oral cavity

Locator R-Tx Abutment*		
	transgingival height	
	1.00	31110-01
	2.00	31110-02
	3.00	31110-03
	4.00	31110-04
	5.00	31110-05
	6.00	31110-06

* All-in-One packaging that includes an abutment, a spacer ring, a cap and an insert for each retention (gray, blue, pink, white).

Accessories for Locator R-Tx abutment		
	unique platform	
	30013	pack of 4 caps and 4 black retainers for impression taking
	30001	grey, without retention, pack of 4 pcs
	30002	blue, low retention, pack of 4 pcs
	30003	pink, medium retention, pack of 4 pcs
	30004	white, high retention, pack of 4 pcs
	30012	black, for the impression taking, pack of 4 pcs
	8514	spacer rings for the prosthesis relining phase, pack of 20 pcs
	30014	aluminium analogs for abutments ø 3.35 mm pack of 4 pcs
	30015	aluminium analogs for abutments ø 4.00 mm, pack of 4 pcs
	30016	aluminium analogs for abutments ø 5.00 mm, pack of 4 pcs
	30017	transfers with low retention, pack of 4 pcs
	30018	spacer rings for relining phase, pack of 4 pcs
	30021	Locator Core Tool for Locator R-Tx
	08008	RA screwdriver for Locator R-Tx, hex 1.25 mm
	08007	digital screwdriver for Locator R-Tx, hex 1.25 mm

Ball attachments and accessories

Ball attachments		
	transgingival height	unique platform
	1.00	A-AS-330-1
	2.00	A-AS-330-2
	4.00	A-AS-330-4

Caps for ball attachments		
	unique platform	
	CAP-TFL-1	polyamide cap
	CONT-CAP-TFL-1	steel container for polyamide cap
	CAP-TIT-1	Gr. 5 Ti cap with cap in two parts, Ti retention spring and plastic
	AN-CAP-TIT-1	spare plastic ring
	MOL1-CAP-TIT-1	spare steel spring, average hardness
	MOL2-CAP-TIT-1	spare steel spring, soft
	CAP-1	gold alloy cap, with 3 plastic rings and a tin spacer for the lab

Components for bars

Overdenture on bars		
	BARC-CAV-TIT	castable bar, l 5.00 cm, h 3.00 mm, thickness 2.20 mm, ovoid-shaped profile with spacer
	CAV-TIT	divisible bar attachment in titanium for oval bars, h 3.00 mm, thickness 2.20 mm
	BARC	castable bar, l 5.00 cm, ø 2.20 mm
	CAV-375	bar attachment in gold alloy for round bars ø 2.20 mm

Analog for ball attachments	
	unique platform
	ANAS

Blank lined paper with horizontal ruling lines.



rev. 07-22



Sweden & Martina S.p.A.

Via Veneto, 10
35020 Due Carrare (PD), Italy
Tel. +39.049.9124300
Fax +39.049.9124290
info@sweden-martina.com
www.sweden-martina.com

Sweden & Martina Ltd

Unit 1b
Amberley Court, Whitworth Road
Crawley, West Sussex, RH11 7XL
Toll free 0800 1123575
info.uk@sweden-martina.com

Sweden & Martina Inc.

Distributor for U.S.
info.us@sweden-martina.com

Sweden & Martina Mediterranea S.L.

España - info.es@sweden-martina.com

Sweden & Martina Lda

Portugal - info.pt@sweden-martina.com

The implants, standard prosthetic components and surgical instruments contained in this catalogue are medical devices and are manufactured by Sweden & Martina S.p.A., except for:

- the Locator Abutment and related components which are medical devices manufactured by Zest Anchors, Inc., 2875 Loker Avenue East, Carlsbad, CA 92010, USA. Locator is a trademark by Zest Anchors, Inc.
- the Dynamic Abutment, Interfase Dynamica which are medical devices manufactured by Talladium España S.L., Avenida Blondel, 54 3º, 25002 Lleida, Spagna. Dynamic Abutment and Interfase Dynamica are trademarks by Talladium.

They conform to the ISO 9001 and ISO 13485 standards and are certified with the CE Mark (Class I) and CE 0476 mark (Class IIA and class IIB) in compliance with Regulation (EU) Medical Devices n. 2017/745.

The Sweden & Martina plant manufactures Medical Devices in compliance with the CGMPs in force in the USA and in other countries worldwide.



Some products may not be regulatory/released for sale in all markets.

All trademarks herein are the property of Sweden & Martina S.p.A. unless otherwise indicated.

This material is intended for laboratories and clinicians and is not intended for patient distribution.

This material is not to be redistributed, duplicated, or disclosed without the express written consent of Sweden & Martina S.p.A.

For additional product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see Sweden & Martina S.p.A. website.

The contents are updated at the time of publication. Check with the company for any subsequent updates.